



User Guide for Coccyx Seat Cushion by AGLLUX

* How to use the seat cushion

- Locate the seat cushion in the middle of a seat or chair as shown in the picture below.



- The cutout for the tailbone should be placing facing the back of your seat or chair.
- We recommend that the coccyx seat cushion be used by people with a weight of no more than 265 pounds (120 kg).

*** Cleaning and washing the seat cushion**

- Locate the zipper on the coccyx seat cushion. Unzip the covering and carefully remove the outer upholstery. It is recommended to wash the upholstery in water up to 104° F (40°C).
- Once clean, dry the upholstery and place it back over the memory foam.
- Do not wash the memory foam in a washing machine! To clean the memory foam, dampen a sponge and carefully remove the dirt. Also, we do not recommend drying the memory foam in a dryer!

*** Customer Support**

- If you have any questions, please contact us prior to contacting Amazon.
- We provide a 1-year warranty on our products. If any of our products that you bought are faulty or turn out not to serve your needs well, please email us. We will provide a suitable replacement or refund! You will not even have to send the products back!

Email: support@agllux.com

Best regards,

Team Support of AGLLUX

Disclaimer

The information provided by AGLLUX is not intended as medical advice its intention is solely informational, and educational. You should consult a health care professional should the need for one be indicated. If you have any concern regarding your health, consult your doctor before use. Use of the seat cushion, or the information provided by AGLLUX should not be a substitute for professional advice, and / or consultation with a licensed health care professional. Do not use the seat cushion to attempt treatment of serious ailments, health conditions, or injuries without the advice of your doctor. No medical claims are warranted or implied by the use of our products. AGLLUX will not be held liable for misuse of the seat cushion in attempted treatment of unknown origins, or undetermined pain. Consult your doctor for treatment of any health condition before beginning any exercise program or initiating use of any physical support product.