U.K. Patent No. 2057016
U.S. Patent No. D379,315 U.K. Patent No. 2057016

**Learning
Resources®**

Your opinion matters! Visit
LearningResources.com
to write a product review or to
find a store near you.

Like
us on


ATENCIÓN: PELIGRO DE ASFIXIA.
Piezas pequeñas. No se recomienda para menores de 3 años.
ATTENTION: RISQUE D'ÉTOUFFEMENT.
Petites pièces. Interdit aux enfants en dessous de 3 ans.
ACHTUNG: ERSTICKUNGSGEFAHR.
Kleine Teile. Nicht geeignet für Kinder unter 3 Jahren.

© Learning Resources, Inc., Vernon Hills, IL, US
Learning Resources Ltd., Bergen Way,
King's Lynn, Norfolk, PE30 2JG, UK
Please retain the package for future reference.
Made in China. LRM2420-GUD
Hecho en China. Conserva el envase para
futuras consultas.
Fabriqué en Chine. Veuillez conserver
l'emballage.
Hergestellt in China. Bitte Verpackung gut
aufbewahren.

Precision School Balance with Weights

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Activity Guide

Guía de actividades • Guide d'activités • Spielvorschläge

 **WARNING:**
CHOKING HAZARD - Small parts.
Not for children under 3 years.

 **CAUTION:**
Federal law requires us to advise that THE BRASS WEIGHTS
in this educational product may contain lead which may be
harmful if swallowed. We stand behind the safety of all of our
products. For further information please contact us at
QA@learningresources.com or 800-222-3909.

Help students explore, develop, and extend their understanding of mathematical and scientific concepts including inequalities, mass, comparisons, the four operations, and order. The Precision School Balance comes with a 10-piece gram set including two 1-gram masses, two 2-gram masses, two 5-gram masses, two 10-gram masses, one 20-gram mass, one 50-gram mass.

Balance Components

- The damping button helps make quick and accurate measurements. The button steadies the swing arm and prevents the pans from rocking during a measurement.
- The calibration slider allows you to zero the balance.
- Marked calibrations demonstrate when the balance is level and zeroed.
- The balance contains a case with a lid to hold the masses. Storing your masses in this compartment will reduce the likelihood of misplacing them.

Use the balance in directed activities such as those listed below or simply place it in a learning center for students to explore. Either way, this balance will expand and enrich the learning opportunities in your classroom.

Exploring the Precision School Balance

Collect a variety of classroom objects such as pencils, erasers, paper clips, scissors, counters, or blocks. Encourage students to place objects on either side of the balance and observe the results. Ask questions that will increase your students' explorations.

You might ask:

- Are there any objects that balance each other?
- Does size always tell which object has the greatest mass?
- Can a small object have greater mass than a large object?

After students have explored with objects on either side of the balance, encourage them to explore measurement with the unit masses provided with the balance. Students can place one object on one side of the balance and use masses to find out if an object has greater, lesser, or close to equal the mass in a unit weight.

Equal Masses

Once students have explored measurements that are greater, lesser, or close to the actual mass, encourage them to find the actual mass. Introduce them to the calibration lines that show when the balance is level. Ask, "Which line should the pointer be on for objects to be balanced?" (Answer: The center line.)

Discuss how to determine an object's mass. Provide examples and allow students to discover the most efficient method. You might prompt their explorations with questions such as:

- Is it best to start with the small unit masses to measure an object?
- Is it best to start with the large unit masses to measure an object?

Ordering Masses

Encourage students to arrange, by mass, the items they have collected. Students should choose one item, determine its mass, and record its mass in a data table. After doing this for all the objects, students should arrange them in order from the least mass to the greatest mass.

Mass Hunt

Play a game of Mass Hunt! Hold up an object and ask students to identify it, and then determine its mass. Challenge them to find another object in the classroom that has a similar mass. Give each student the opportunity to choose and reject items until they agree on one item that they think has the mass that most closely equals that of the object you selected. Discuss students' choices and then determine the masses of the most popular choices to see which one is the best answer.

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