smart
weigh

Digital Body Fat Scale Manual

Copyright © 2016 Smart Weigh™

Index

English 4

Español 15

Français 21

Italia 27

Deutsch 33

Getting Started

Our team at Smart Weigh thanks you for purchasing the Smart Memory Scale. Our goal is to bring you precision measuring devices that will ensure many years of quality and reliability. Please read through the safety and operating guidelines for optimal performance. Keep this manual and warranty for future reference.

- Always weigh on a flat and leveled surface, free from vibrations and drafts. The corner of a room is usually the most stable.
- Avoid exposure to extreme heat or cold. Scales perform best at normal room temperature. If temperatures have changed dramatically, calibration and adjustment may be necessary.
- Do not exceed the scale's maximum capacity. Overloading the scale can permanently damage it!
- Store the scale in a clean, dry location. Dust, dirt, or moisture can accumulate in the weighing sensors causing inaccuracies or malfunction.
- Avoid using the scale in close proximity to radio wave emitting devices such as computers, cash registers, or cordless phones.
- Do not drop items onto the weighing platform. Gently apply all items to be weighed.
- Avoid dropping the scale. The warranty does not cover damage due to rough treatment or overload.
- Do not disassemble the scale. This product contains no user serviceable parts.
- If you are having trouble with your scale, check the battery. This simple step can remedy most scale issues.
- Only operate the scale on a stable, flat surface. Carpeted or uneven surfaces may affect the accuracy.

Overview

- The Smart Weigh Digital Body Fat Scale uses Bio Impedance Analysis (BIA) technology which sends a low-level electrical current through the body to calculate body fat, body water, muscle mass, and bone mass. The electrical current is mild and may not be felt. Contact is made through the scale platform.
- The information provided by this device is not meant to treat, cure, or prevent any disease or illness. This device should not be used by anyone who is acutely or chronically ill, suffering from a disease, or taking medication that affects your water levels. The accuracy for such patients has not been verified. Specific medical advice should be obtained from a physician.

- This scale is not intended for use by pregnant women or children under the age of 10. It should not be used by anyone with an implanted medical device (pacemaker, metal plates or screws, or contraceptive devices). If in doubt, contact your physician.
- **Note: Bare feet are required for the scale to function properly.**

Battery Installation

- Before using the scale for the first time, check the battery compartment and remove any plastic insulation strips or packaging that may be in place to prevent battery drainage during shipping and storage.
- In the event that the LCD displays "Lo", replace with 2 x CR2032 lithium batteries. Be sure to properly align the positive and negative contacts.

Changing Modes

1. Once the scale is on, press the button on the bottom of the scale to select the desired unit measurement.
2. Set the scale for kilograms (kg), stones (st/lb), or pounds (lb).
3. The scale needs to be initialized before the first use and after battery replacement. To initialize, press firmly on the scale platform. The display will show a running zeros pattern and automatically turn off.

Personal Profile Setup

- To use the body fat analyzer you must first set up a profile with your personal information so that your body composition can be calculated.
- You must be barefoot when taking a reading. The body fat measurements can only be calculated when your skin comes in direct contact with the scale platform.

Please note: During user setup, if the scale is inactive for 20 seconds, the scale will divert to general weighing mode and user data will not be saved.

1. The scale can store up to 8 user profiles. They can be used to store information such as gender, height, age, and weight for each user. This scale provides measurements of body fat, body water, bone mass, and muscle mass which are calculated according to the users personal data.
Note: If you only desire to use the scale for general body weighing you do not need to set up a profile.

2. To begin setting up your profile, press the SET button. The scale will enter the data input mode.
3. Press the **▲** or **▼** arrow to select a user number from the list of available users. All your information will be stored under this user ID number. Press SET to confirm.
4. Press the **▲** or **▼** arrow to select the correct gender and athlete mode. Choose between male, female, male athlete, and female athlete icons. Press SET to confirm.
5. Press the **▲** or **▼** arrow to increase or decrease the height to the correct value. Press SET to confirm.
6. Press the **▲** or **▼** arrow to increase or decrease the age to the correct value. Press SET to confirm.
7. After all settings are entered, scale will divert to weighing mode. "0.0" will be displayed on the screen. Step on the scale immediately to take the first weight reading.

Note: To complete user setup, the initial weight reading must be taken. If user does not step on the scale or scale is inactive for 5 seconds, the scale will automatically turn off and no information will be saved.

8. Weight will stabilize on the screen and the scale will proceed to calculate other body measurements. While analysis is in progress the screen will display a running zeros pattern **0 0 0**.
9. Next, your readings of body fat, body water, muscle mass, and bone mass will show on the screen.
10. To create additional users, repeat the steps above.

Scale Operation

General Weighing

The scale can operate as a conventional weight-reading scale without any special programming steps required. To use the scale in general weighing mode only:

1. Remove any plastic insulation strips or packaging from the battery compartment.
2. Once the scale is initialized, position the scale on a hard flat surface as carpeted or uneven surfaces may affect the accuracy.
3. Step onto the scale platform and stand still. After a few seconds, you will see the weight flashing to indicate a weight reading.
4. Your weight will display for 5 seconds and the scale will automatically shut off.

Body Fat Analyzer

If the scale recognizes the user and the user is barefoot on the platform, the scale will automatically display the user's body fat readings. To set up a profile and become a recognized user please see the steps outlined in the "Personal Profile Setup" section.

1. Step onto the scale platform with bare feet and remain still while the scale computes your weight. Your weight will then appear.
2. Based on your previous weight, the scale will recognize you as a user and will calculate your body fat readings according to the data saved in your profile. While weight is stored and analysis is in progress the screen will display a running zeros pattern **0 0 0**.
3. When the scale successfully recalls your user profile, it will display the body fat, water, muscle, and bone results in order.
4. In the event that two users are within the same weight range the scale will display the User ID numbers of each of those users on the screen. The user must then select the correct user number to identify himself by pressing the **▲** button to select P1 or the **▼** to select P2.
5. If three or more users have similar weights, the scale will scroll two numbers per screen. To select the correct user profile, wait until the correct screen appears, then press the **▼** to select the number on the left or the **▲** to select the number on the right.

To select your User ID manually follow the steps below:

1. Press SET to turn the scale on.
2. Press the **▲** or **▼** button to select your personal user ID. Your personal data will appear on the screen.
3. Wait until 0.0 appears on the screen and then step on the scale with bare feet. Stand still while the weight is measured.
4. Your weight will show for two seconds.
5. Continue to stand still on the scale while your other readings are calculated. The screen will show a running zeros pattern **0 0 0** while the scale computes.
6. The body fat, water, muscle, and bone readings appear next.
7. The results will show three times and then the scale will turn off.

Troubleshooting

Problem	Solution
Scale will not turn on	Check the batteries
Displayed weigh fluctuates randomly	Place scale on flat surface
Displayed weight is inaccurate	Recalibrate the scale by stepping on the scale and stepping off.
	Clothes can cause inaccuracies in weight readings.
	Food and drink consumption prior to weighing can cause inaccuracies.
	Scale is bumped while in use. Move scale to stable location.
	Place your feet in the same position at every weight reading
Scale will not display body composition measurements.	Be sure to set up your personal profile before taking a reading
	Remove shoes and socks from feet before stepping on the scale. Be sure that your skin comes in direct contact with the surface of the scale.
	Dampen the bottom of your feet before stepping on the scale.

Error Codes

1. "----" – Maximum Capacity Exceeded. Remove the extra weight immediately as permanent damage can occur to the load cell.
2. "Lo" – Batteries are low. Please replace with 4 new AAA batteries.
3. If the weight is all that appears after the Bio Impedance Analysis (BIA) was attempted, or you get Instability Error, Contact Error, or Body Fat Range Error, repeat the measurement again while maintaining maximum contact between your feet and metal sensors.

Please note: If the scale has been moved or has been inactive for two days, step onto the scale and step right off. This process ensures the scale is set and ready to use. At other times step right onto the scale for weight readings.

Care and Maintenance

- The product is intended for home / consumer use only; it is not intended for professional use in hospitals or medical facilities.
- The body fat percentage of a pregnant woman could be inaccurate and misleading due to the changes in the body composition. Pregnant women should consult their doctor before use.
- This scale uses BIA (Bioelectrical Impedance Analysis) to determine body fat percentage. BIA sends a harmless signal through your body. Do not use this product if you have a pacemaker or any other internal medical device. When in doubt, contact your physician.
- Clean with a lightly dampened cloth. Do not use solvents or immerse the product in water. Thoroughly dry the scale after cleaning as moisture can damage the scale's sensors and electronics.
- Do not use the scale on carpet or uneven surfaces and do not drop or jump on it. Treat it with care to ensure the best performance.
- Do not store the scale where you store cleaning chemicals. The vapors may affect the electronic components of your scale. Do not store the scale on its side.
- Store your scale in a clean, dry location at room temperature. Dust, dirt, and moisture from humidity can accumulate on the weighing sensors causing inaccuracy or malfunction.
- To prevent battery drainage, do not store anything on the scale. Do not dispose of batteries in fire. Batteries may explode or leak.
- Do not compare weight readings from one scale to another as some differences will exist due to manufacturing tolerances.
- It is recommended that you measure the same time of the day, preferably early evening before a meal, for the most consistent results.
- Do not disassemble this products other than replacing batteries. It contains no user serviceable parts.
- If you weigh yourself twice and get two different readings, your weight lies between the two readings.

Body Mass Index

Your BMI is calculated based on your height and weight. As your body fat, water, muscle and bone results appear, your BMI appears as a symbol in the bar at the bottom of the display. The symbols are:

- Indicates Underweight

+ Indicates Overweight

0 Indicates Normal Weight

++ Indicates Obesity

Maintaining a proper weight is essential to healthy living. Obesity has been linked to a higher risk of heart disease, diabetes, and some forms of cancer. The Body Mass Index (BMI) is calculated based on a person's height and weight, and is the most widely accepted measurement of weight according to health professionals. Use the chart below to determine your own BMI by matching your height in the left column with your weight in the center. Anyone, male or female, with a BMI of 25–29 is considered overweight and a body mass index of 30 is considered obese.

BMI	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35
Height	Weight in Pounds																
4'10"	91	96	100	105	110	115	119	124	129	134	138	143	148	153	158	162	167
4'11"	94	99	104	109	114	119	124	128	133	138	143	148	153	158	163	168	173
5'	97	102	107	112	118	123	128	133	138	143	148	153	158	164	169	174	179
5'1"	100	106	111	116	122	127	132	137	143	148	153	158	164	169	174	180	185
5'2"	104	109	115	120	126	131	136	142	147	153	158	164	169	175	180	186	191
5'3"	107	113	118	124	130	135	141	146	152	158	163	169	175	180	186	191	197
5'4"	110	116	122	128	134	140	145	151	157	163	169	174	180	186	192	197	204
5'5"	114	120	126	132	138	144	150	156	162	168	174	180	186	192	198	204	210
5'6"	118	124	130	136	142	148	155	161	167	173	179	186	192	198	204	210	216
5'7"	121	127	134	140	146	153	159	166	172	178	185	191	198	204	211	217	223
5'8"	125	131	138	144	151	158	164	171	177	184	190	197	203	210	216	223	230
5'9"	128	135	142	149	155	162	169	176	182	189	196	203	209	216	223	230	236
5'10"	132	139	146	153	160	167	174	181	188	195	202	209	219	222	229	236	243
5'11"	136	143	150	157	165	172	179	186	193	200	208	215	222	229	236	243	250
6'	140	147	154	162	169	177	184	191	199	206	213	221	228	235	242	250	258
6'1"	144	151	159	166	174	182	189	197	204	212	219	227	235	242	250	257	265
6'2"	148	155	163	171	179	186	194	202	210	218	225	233	241	249	256	264	272
6'3"	152	160	168	176	184	192	200	208	216	224	232	240	248	256	264	272	279
	Healthy Weight						Overweight						Obese				

Body Fat

The human body consists of a percentage of fat. Body fat is an important component of a healthy and functioning body because it protects vital organs, helps regulate body temperature, stores vitamins, and helps the body sustain itself when food is scarce. However, too much body fat or too little body fat can cause damage to your health. It is difficult to gauge how much body fat we have in our bodies simply by looking in the mirror. That is why it is important to measure and monitor your body fat percentage. Body fat percentage gives you a better measurement of fitness than weight alone; the composition of your weight loss could mean you are losing muscle mass rather than fat.

Ideal body fat content is not the same for all people. Age, sex, and hereditary are variables in this measurement. Refer to the table below as a general guide. Consult your physician to determine what is ideal for you.

Standard for Men						Standard for Women					
Rating	Age					Rating	Age				
	20-29	30-39	40-49	50-59	60+		20-29	30-39	40-49	50-59	60+
Low	<13	<14	<16	<17	<18	Low	<13	<14	<16	<17	<18
Normal	14-20	15-21	17-23	18-24	19-25	Normal	14-20	15-21	17-23	18-24	19-25
Moderately High	21-23	22-24	24-26	25-27	26-28	Moderately High	21-23	22-24	24-26	25-27	26-28
High	>23	>24	>26	>27	>28	High	>23	>24	>26	>27	>28

It is best to measure body fat and weight daily at the same time. For the best accuracy use this scale without clothing. Feet must be bare and clean and for best readings they should be slightly damp.

Body Water

Body water is an important component of body weight. It represents over half of your total weight and almost two thirds of your lean body mass (predominantly muscle).

Water is critical to healthy function in the body. It performs a number of important roles in the body. All the cells in the body can only function properly if they have enough water. Water also plays a vital role in regulating the body's temperature, specifically through perspiration. The combination of your weight and fat measurements may appear to be normal, but your body hydration level could be insufficient for healthy living.

This scale measures your weight, the electrical impedance between your feet, and then combines those readings with personal information you entered (height, age, gender). It uses this information to calculate your body fat percentage. Please note that your level of hydration (water content) varies throughout the day and thus affects your impedance measurements. If you are dehydrated, there is a greater chance of fluctuation in the body fat percentage displayed by this scale than when you are fully hydrated.

For best results, choose a consistent time to use this analyzer, such as in the morning. By establishing the pattern of your body's fluctuations over a period of one month, you can determine your average level. Women will have more fluctuations in water content than men and need to monitor their readings over a longer period to establish a personal average. Wait several hours before taking a reading when drinking coffee or alcohol, taking diuretic medications, or exercising. These all affect the level of hydration and the accuracy of the scale.

	% Body Fat Range	Normal % TBW Range
Men	14 to 14 %	70 to 63 %
	15 to 21 %	63 to 57 %
	22 to 24 %	57 to 55 %
	25 and over	55 to 37 %
Women	4 to 20 %	70 to 58 %
	21 to 29 %	58 to 52 %
	30 to 32 %	52 to 49 %
	33 and over	49 to 37 %

Body water measurements results are influenced by the proportion of body fat and muscle. If the proportion of body fat is high, or the proportion of muscle is low, then the body water results will tend to be low.

Measurements such as body weight, body fat, and body water are tools for you to use as part of your healthy lifestyle. Since short term fluctuations can be normal, we suggest that you chart your progress over time, rather than focus on a single day's reading.

Muscle Mass

An inactive lifestyle can cause a loss of lean body mass while body fat increases. There are many benefits in gaining muscle mass such as: reversing the decline in strength, bone density, and muscle mass with aging, and maintenance of flexible joints.

Proper calorie intake and nutrition are vital in maintaining a healthy body. Consult your physician before undertaking any exercise or diet program.

Bone Mass

Regular exercise and a balanced diet can help maintain healthy bones. Like muscle, bone is a living tissue that can respond to exercise by becoming stronger. For most people, bone mass peaks in their thirties. Then people begin to lose bone. A regular exercise regimen will help prevent that loss. Calcium and vitamin D contribute to healthy bones.

Bone mass readings given by this product are an estimation of the amount of bone in your body. Individuals with osteoporosis or low bone densities may not get accurate estimations. Please consult your doctor with any concerns.

SBS500 at a Glance

Durable quality ensures years of reliable use

Measures body weight along with BMI, body fat, body water, bone mass, and muscle mass

Accurate, consistent and fast weight readings displayed on large LCD screen

Modern design combines reflective black tempered glass surface with stainless steel accents.

Step-on technology triggers power upon contact

Stores up to 8 user profiles with automatic recognition

Power saving automatic shut off

Capacity	440 lbs / 200 kg / 31 st: 6 lbs
Readability	0.2 lb / 0.1 kg / 1lb(0.2lb)
Dimensions	12.6" x 12.6" x 0.75"
Modes	kg / lb / st;lb cm / in
Body Fat Graduation	0.1%
Height Range	3'3.5" - 7'2.5" (100 - 220 cm)
Age Range	10+
Measurements	weight, bmi, body fat, body water, muscle mass, water mass
Users	8
Power	2 x CR2032 lithium batteries
Weight	4.4 lbs / 2 kg

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
1. Press SET to turn the scale on.
2. Press  or  button to select a desired memory number. Press the SET button to confirm. The gender/athlete mode blinks.
3. Press  or  button to switch between the male, female, male athlete, and female athlete icons. Press SET to confirm your gender / athlete selection. The height digits blink.
4. Press  or  to increase or decrease the height value. Press the SET button to confirm the displayed height. The age digits blink.
5. Press  or  button to increase or decrease the age value. Press the SET button to confirm the displayed age.

When the scale shows 0.0, immediately step onto the scale for a measurement and remain standing. Otherwise, the scale will turn off automatically. The scale will first display your weight and the data will be saved. Next your body fat, water, muscle, and bone will show on the screen.

If there is no operation after five seconds and you have not completed the set-up, the scale will turn to 0.0 and the process will have to be done again.

Use the above procedure to set any additional users.

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Warranty

The Smart Weigh Scale is warranted to the original owner for two years from the date of purchase for defects in material and workmanship under normal use and service. If your scale proves to be defective within one year of the purchase due to a manufacturer's defect, then follow the terms outlined below:

- Pack your scale carefully to avoid damages in transit. This warranty will not apply to damages that occurred in the return transit.
- This warranty will not cover any product that has been used for commercial purposes, rental use, or use that is not intended for this scale.
- Smart Weigh will not be liable for any damages, losses, or inconveniences caused by user negligence, abuse, or use that is not in accordance with this product manual.
- Smart Weigh is not responsible for any incidental or consequential damages of any nature resulting from use of this product. Any liability shall not exceed the purchase price of this product.

Smart Weigh Customer Support

Email: support@smartweighscales.com

Website: www.smartweighscales.com

Phone: 855.797.2253