## Quick Start

### Plan your Veet® InfiniSilk™ Pro™ hair removal treatments for best results

Treatments 1-4 – plan two weeks apart

Treatments 5-7 – plan four weeks apart

Treatments 8 + – treat as needed, until desired results are achieved.

### Set up Veet® InfiniSilk™ Pro!

1. Remove Veet® InfiniSilk™ Pro™ Pulse Light Technology Device and other components from box.
2. Verify that the DISPOSABLE LAMP CARTRIDGE is inserted properly at the tip of the device.
3. Plug the ADAPTOR cord into the Veet® InfiniSilk™ Pro™ HAIR REMOVAL UNIT socket.
4. Plug the ADAPTOR into an electrical outlet. The STAND BY INDICATOR LIGHT will turn on. Your Veet® InfiniSilk™ Pro™ is now ready to start.

### Treat with Veet® InfiniSilk™ Pro safely and properly!

5. Choose areas to treat that have not had recent exposure to the sun.
6. The skin of the treatment areas should be shaved, clean, dry and free of any powders, antiperspirants or deodorants. Do not wax, pluck or tweeze the hair out.  
Shaving is an important step to get the results you want!
7. Press the ON/STBY SWITCH. The POWER ON INDICATOR LIGHT (orange) will turn on and a fan sound (similar to the sound of a hairdryer) will start.
8. Approximately 3.5 seconds after pressing the ON/STBY SWITCH, the READY INDICATOR LIGHT (green) will turn on. The device is then ready for you to trigger the first pulse at the lowest energy level.
9. If this is your first treatment - Choose the energy level.  
InfiniSilk™ Pro offers 5 levels of energy from 1 (the lowest) to 5 (the highest).  
To choose the right energy level for your first treatment, test each treatment area separately:

- Place the TREATMENT SURFACE on the skin and press the trigger to deliver one pulse.
- If you did not experience abnormal discomfort press the + button once and deliver one pulse on a different spot at energy level 2.
- Continue testing at increased energy levels until you have reached the highest level that you feel comfortable with. If within one hour you experience no adverse effects, you can start your full treatment at that energy level.
- Repeat this test for each body part you intend to treat.

10. Set the energy level. If the energy level should be higher than the lowest, press the ENERGY LEVEL SETTING BUTTONS using the “-” or “+” to respectively decrease or increase the energy level, until the desired energy level is set and represented by 3 energy level indicator lights (For details see “Energy Level” chapter inside this user manual).
11. Apply the TREATMENT SURFACE to the skin, making sure the skin is spread evenly and smoothly. As soon as the TREATMENT SURFACE is in full contact with the skin, the READY INDICATOR LIGHT will start blinking.
12. Press the PULSE BUTTON. The device will first determine the color of your skin. If the color of the skin is light enough for safe application, the device will flash a pulse of light onto your skin, and the READY INDICATOR LIGHT will turn off.  
You will see a bright flash of light and simultaneously hear a subtle pop sound, which is a normal noise for the device. You will feel a mild sensation of warmth and tingling.  
Veet® InfiniSilk™ Pro™ will immediately recharge for the next pulse. After 3.5 seconds the READY INDICATOR LIGHT will turn on again.
13. Remove the TREATMENT SURFACE from the treated area of skin.
14. Move the TREATMENT SURFACE to another area of skin. Use the pressure marks the device just made on your skin to guide you for proper positioning of the next pulse, avoiding both gaps and overlaps between pulses.

Try to avoid overlapping pulses!

**Warning:** Do not treat the same area of skin more than once per hair removal session!  
Treating the same area of skin more than once per session increases the likelihood of adverse effects.

15. Repeat the process until the entire treatment area is covered.

If your skin blisters or burns, STOP USE IMMEDIATELY!

Copyright 2014 © RB. All rights reserved.  
Print date: February 2014

Distributed by : Reckitt Benckiser LLC

Parsippany, NJ 07054-0224

For questions and concerns, call 1-866-252-5327

www.veet.us

“VEET, VEET logo and INFINI’SILK are trademarks of the Reckitt Benckiser Group of Companies.

Veet® and Veet® logo are trademarks of Home Skinovations Ltd, used under license.

## Table of Contents

1. Before You Start .....	4
1.1 Veet® Infini’Silk™ Pro™ device description .....	4
1.2. Package contents .....	5
1.3. The Lamp Cartridge .....	5
2. Veet® Infini’Silk™ Pro™ Intended Use .....	5
3. Safety with Veet® Infini’Silk™ Pro™ .....	5
4. Contraindications .....	6
5. Do not use Veet® Infini’Silk™ Pro™ if: .....	9
6. Avoid using Veet® Infini’Silk™ Pro if any of the following applies: .....	9
7. Possible Side Effects using Veet® Infini’Silk™ Pro™ .....	10
8. Long Term Hair Removal the Veet® Infini’Silk™ Pro™ Way .....	11
9. Get Started! .....	12
9.1. What to Expect when treating with Veet® Infini’Silk™ Pro™? .....	16
9.2. After treating with Veet® Infini’Silk™ Pro™ .....	16
10. Maintenance of Veet® Infini’Silk™ Pro™ .....	17
10.1. Cleaning Veet® Infini’Silk™ Pro™ .....	17
10.2. Replacing the LAMP CARTRIDGE .....	17
11. Troubleshooting .....	18
11.1. “My Veet® Infini’Silk™ Pro™ does not start.” .....	18
11.2. “A light pulse is not emitted when I press the PULSE SWITCH” .....	18
12. Consumer Relations .....	19
13. Specifications .....	19
14. Labeling .....	59

Before using Veet® InfiniSilk™ Pro™ for the first time, please read this User Manual in its entirety. Pay particular attention to sections on device use procedures, device operation, and after-use procedures. We recommend you re-familiarize yourself with this User Manual before each use of Veet® InfiniSilk™ Pro™. Veet® InfiniSilk™ Pro™ is a powerful electrical device. As such, it should be used with special attention to safety.

Please read all warnings and safety precautions before use, and strictly follow them when using Veet® InfiniSilk™ Pro™

## I. Before You Start

Veet® InfiniSilk™ Pro is a light-based device for permanent hair reduction results, designed for home-use.

### 1.1 Veet® InfiniSilk™ Pro™ device description

Your Veet® InfiniSilk™ Pro™ consists of a HAIR REMOVAL DEVICE with a CONTROL PANEL and a REPLACEMENT LAMP CARTRIDGE.

On the Veet® InfiniSilk™ Pro™ HAIR REMOVAL UNIT HANDLE you can find the PULSE BUTTON. The DISPOSABLE LAMP CARTRIDGE is located at the tip of the HAIR REMOVAL UNIT.

- ① Hair Removal device
- ② Control panel
- ③ Pulse button
- ④ lamp cartridge
- ⑤ Skin Color Sensor
- ⑥ Treatment Surface
- ⑦ Control Panel
- ⑧ ON/STBY switch
- ⑨ Stand by indicator light
- ⑩ Power on indicator light
- ⑪ Energy level setting buttons
- ⑫ Energy level indicator lights
- ⑬ Replace cartridge indicator light
- ⑭ System warning indicator light
- ⑮ Skin tone warning indicator light
- ⑯ Ready indicator light

### 1.2. Package contents

Upon opening the Veet® InfiniSilk™ Pro™ package, you will find the following parts:

- Veet® InfiniSilk™ Pro™ Hair Removal device
- A wall mount adaptor
- This User's Manual

### 1.3. The Disposable Lamp Cartridge

The Veet® InfiniSilk™ Pro™ all-purpose LAMP CARTRIDGE can fire a limited and pre-set number of light pulses (regardless of the energy level of these pulses).

The type of LAMP CARTRIDGES included with your system, and the number of light pulses in them, is specified on the package of your Veet® InfiniSilk™ Pro™.

## 2. Veet® InfiniSilk™ Pro™ Intended Use

Veet® InfiniSilk™ Pro™ is intended for removal of unwanted hair. Veet® InfiniSilk™ Pro™ is also intended to effect long term, or permanent hair reduction. Permanent hair reduction is defined as a long-term stable reduction in the number of hairs re-growing after a treatment regime.

Veet® InfiniSilk™ Pro™ may be used to remove unwanted body hair.

Ideal body areas for Veet® InfiniSilk™ Pro™ use include the underarms, bikini line, stomach, arms, legs, back shoulders and chest.

## 3. Veet® InfiniSilk™ Pro - Safe to Use (as directed)

With Veet® InfiniSilk™ Pro™ Safety Comes First.

### ■ Veet® InfiniSilk™ Pro - Sensor for the skin

Light based hair removal is not suitable for naturally darker skin tones or tanned skin. Veet® InfiniSilk™ Pro™ comes with a built-in SKIN COLOR SENSOR that is designed to measure the complexion of the applied surface and enable application only on suitable skin complexions. This safety feature will not let you treat where your skin is too dark or too tanned.

Furthermore, the 1.3X3cm<sup>2</sup> TREATMENT SURFACE through which pulses of light are delivered is recessed inside the LAMP CARTRIDGE. This enables Veet® InfiniSilk™ Pro™ to protect your skin by avoiding direct contact between the TREATMENT SURFACE and the skin.

#### ■ Veet® InfiniSilk™ Pro - No need for eyewear

Veet® InfiniSilk™ Pro™ has a built-in safety feature for eye protection. It has been designed so that a light pulse can not be emitted when the TREATMENT SURFACE is facing open air. The safety switch enables treatment only when the TREATMENT SURFACE is in contact with the tissue.

## 4. Contraindications

### Important Safety Information – Read Before Use!

Veet® InfiniSilk™ Pro™ is not designed for everyone. Please read and consider the information in the following section before use.

For further information and personalized advice you may also visit [www.veet.us](http://www.veet.us).

If you are unsure if Veet® InfiniSilk™ Pro™ is safe for you to use, please ask your doctor or dermatologist!

Warning	Possible results	Safety Tip
<p><b>Do not use Veet® InfiniSilk™ Pro on naturally dark skin complexion.</b></p> <p>Not Safe!</p> 	<p>Treating dark skin with Veet® InfiniSilk™ Pro can result in adverse effects such as redness and skin discomfort</p>	<p>A SKIN COLOR SENSOR is embedded in Veet® InfiniSilk™ Pro™ to measure the treated skin complexion at the beginning of each session and occasionally during the session. The SKIN COLOR SENSOR ensures that pulses will only be emitted on suitable skin tones.</p>

Warning	Possible results	Safety Tip
<p><b>Do not use Veet® InfiniSilk™ Pro on tanned skin or after recent sun exposure (including artificial tanning machines). Do not expose treated areas to the sun.</b></p>	<p>Treating with Veet® InfiniSilk™ Pro before or after any sun exposure can result in adverse effects such as redness and skin discomfort. Tanned skin particularly following sun exposure, contains large quantities of the pigment Melanin. all skin types and complexions, including those which don't seem to tan quickly. The presence of large quantities of Melanin exposes the skin to higher risk when using any method of light based hair removal.</p>	<p>Avoid exposure to the sun for 4 weeks before your Veet® InfiniSilk™ Pro™ treatment! Sun exposure includes unprotected exposure to direct sunlight of over 15 minutes constantly, or unprotected exposure to diffused sunlight of over 1 hour constantly. To protect recently treated skin when exposed to sunlight, be sure to thoroughly apply sunscreen SPF 30 or higher, for 2 weeks after each hair removal session.</p>
<p><b>Never use Veet® InfiniSilk™ Pro around or near the eyes.</b></p>	<p>Veet® InfiniSilk™ Pro is not recommended for facial treatment</p>	<p>Veet® InfiniSilk™ Pro™ has a built-in safety feature for eye protection. It has been designed so that a light pulse can not be emitted when the TREATMENT SURFACE is facing open air. The safety switch enables treatment only when the TREATMENT SURFACE is in contact with the tissue.</p>
<p><b>Do not treat on tattoos or permanent makeup, dark brown or black spots, (such as large freckles, birth marks, moles or warts), nipples, genitals or lips.</b></p>	<p>Treating any area with the Veet® InfiniSilk™ Pro that is dark in colour or has more pigment can result in adverse effects such as burns, blisters, and skin color changes (hyper- or hypopigmentation).</p>	<p>Cover area with material that will not absorb the light such as a white cloth or white medical tape.</p>
<p><b>Do not treat on active eczema, psoriasis, lesions, open wounds or infections (cold sores), abnormal skin conditions caused by systemic or metabolic diseases (diabetes for example)</b></p>	<p>Light based treatments can cause adverse reactions on already sensitive areas.</p>	<p>Wait for the effected area to heal before using Veet® InfiniSilk™ Pro™.</p>

Warning	Possible results	Safety Tip
A treatment area that has a history of herpes outbreaks.	Light based treatments can cause adverse reactions on already sensitive areas.	Consult your physician to received preventative treatment before using Veet® Infini'Silk™ Pro™.
Do not treat with Veet® Flash & go if you are pregnant or nursing.	Hormonal changes during pregnancy could increase skin sensitivity and the risk of adverse effects.	
Do not treat on areas where you may want hair to grow.	Results are permanent.	
Keep Veet® Infini'Silk™ Pro away from water! Do not place or store Veet® Infini'Silk™ Pro™ where it can fall, be pushed or placed into a tub, sink or any other vessel containing water.	This may cause severe electrocution.	Do not use Veet® Infini'Silk™ Pro™: while bathing or if it becomes damp or wet. Unplug immediately if it has fallen into water.
Do not attempt to open or repair your Veet® Infini'Silk™ Pro™ device.	Opening Veet® Infini'Silk™ Pro™ may expose you to dangerous electrical components and to pulsed light energy, either of which may cause serious bodily damage and/or permanent eye injury. It may also damage your device and will void your warranty.	Please contact RB Consumer Relations if you have a broken or damaged device in need of repair at 1-866-252-5327 or visit <a href="http://www.veet.us">www.veet.us</a> for questions.

## 5. Do not use Veet® Infini'Silk™ Pro if:

If the LAMP CARTRIDGE is damaged, or if its TREATMENT SURFACE is cracked, chipped or missing.		Stop use immediately and contact RB Consumer Relations.
the attachments or accessories are not recommended by RB.		Stop use immediately and contact RB Consumer Relations.
it has a damaged cord or plug.	Keep the power cord away from heated surfaces.	Stop use immediately and contact RB Consumer Relations.
you see or smell smoke when it is in use.		Stop use immediately and contact RB Consumer Relations.
it is not working properly or if it appears damaged.		Stop use immediately and contact RB Consumer Relations.
its SKIN COLOR SENSOR is cracked, or broken.		Stop use immediately and contact RB Consumer Relations.
the outer shell is cracked or is coming apart.		Stop use immediately and contact RB Consumer Relations.

## 6. Avoid using Veet® Infini'Silk™ Pro if any of the following applies:

History of keloidal scar formation.	Consult with your physician prior to using the Veet® Infini'Silk™ Pro.
If you suffer from epilepsy.	
If you have a disease related to photosensitivity, such as porphyria, polymorphic light eruption, solar urticaria, lupus, etc.	

If you have a history of skin cancer or areas of potential skin malignancies.	Consult with your physician prior to using the Veet® InfiniSilk™ Pro.
If you have received radiation therapy or chemotherapy treatments within the past 3 months.	
If you have any other condition which in your physician's opinion would make it unsafe for you to be treated.	
If your skin was treated with a physician-prescribed drug within the past 6 months, please consult with your physician.	

## 7. Possible side effects using Veet® Infini'Silk™ Pro™

When used according to the instructions, side effects and complications associated with use of Veet® Infini'Silk™ Pro™ are uncommon. However every cosmetic procedure, including those designed for home use, involves some degree of risk. Therefore it is important that you understand and accept the risks and complications that can occur with pulsed light hair removal systems designed for home use.

Adverse Reaction	Degree of Risk
<p><b>Minor Skin Discomfort –</b> Although home pulsed light hair removal is generally very well-tolerated, most users do feel some mild discomfort during use, usually described as being a mild stinging sensation on the treated skin areas. The stinging sensation usually lasts during the time of the application itself or for a few minutes thereafter. Anything beyond this minor discomfort is abnormal and means that either you should not continue to use Veet® Infini'Silk™ Pro™ because you are unable to tolerate the hair removal application, or that the energy level setting is too high.</p>	Minor

Adverse Reaction	Degree of Risk
<b>Skin Redness –</b> Your skin may become red right after using Veet® Infini'Silk™ Pro™ or within 24 hours of using Veet® Infini'Silk™ Pro™. Redness generally clears up within 24 hours. See your doctor if redness does not go away within 2 to 3 days.	Minor
<b>Increased Sensitivity to Skin –</b> The skin of the treated area is more sensitive so you may encounter dryness or flaking of the skin.	Minor
<b>Pigment Changes –</b> Veet® Infini'Silk™ Pro™ targets the hair shaft, in particular the pigmented cells in the hair follicle and the hair follicle itself. Nevertheless there is risk of temporary hyperpigmentation (increased pigment or brown discoloration) or hypopigmentation (whitening) to the surrounding skin. This risk of changes in skin pigmentation is higher for people with darker skin tones. Usually discoloration or changes to skin pigment are temporary and permanent hyperpigmentation or hypopigmentation rarely occur.	Rare
<b>Excessive Redness and Swelling –</b> In rare cases treated skin may become very red and swollen. This is more common in sensitive areas of the body. The redness and swelling should subside within 2 to 7 days and should be treated with frequent applications of ice. Gentle cleansing is OK, but one should avoid exposure to sun.	Rare

## 8. Permanent Hair Reduction with Veet® Infini'Silk™ Pro

A typical full hair growth cycle may take 18-24 months. During this time multiple Infini'Silk™ Pro™ sessions may be required in order to achieve permanent results  
The efficiency of hair removal varies from person to person according to body area, hair color, and how Veet® Infini'Silk™ Pro is used.

**Typical Veet® Infini'Silk™ Pro™ hair removal plan during a full hair growth cycle:**

Treatments 1-4 – plan two weeks apart

Treatments 5-7 – plan four weeks apart

Treatments 8 + – treat as needed, until desired results are achieved.

## Expected Treatment Outcome for Hair Removal

Fitzpatrick skin type	Hair Color	Light Fluence [J/cm <sup>2</sup> ]	Anatomic area of the body <sup>2</sup>	Average Number of Treatments <sup>1</sup>	Effectiveness (% of hair removal/regrowth) <sup>1</sup>
I-VI	Light brown to brown	3-4.5	Lower or upper legs	Session 10-12, approx. four weeks apart.	60% Hair reduction
			Arms		50% Hair reduction
			Bikini Line		60% Hair reduction
			Axilla		55% Hair reduction
I-VI	Brown to black	3-4.5	Lower or upper legs	Session 8-10, approx. four weeks apart.	70% Hair reduction
			Arms		60% Hair reduction
			Bikini Line		70% Hair reduction
			Axilla		65% Hair reduction

- Individual response depends on hair type as well as biological factors that may affect hair growth patterns. Some users may respond quicker or slower than the average number of treatments.
- One cannot expect permanent epilation in a single or even double treatment sessions. The duration of the resting period for hair follicles also depends on the body area.

## 9. Get started!

- Remove Veet® InfiniSilk™ Pro™ HAIR REMOVAL device and other components from box.
- Verify that the LAMP CARTRIDGE is inserted properly at the tip of the device and that it is not damaged.
- Plug the ADAPTOR cord into the Veet® InfiniSilk™ Pro™ HAIR REMOVAL UNIT socket.
- Plug the ADAPTOR into an electrical outlet. The STAND BY INDICATOR LIGHT will turn on. Your Veet® InfiniSilk™ Pro™ is now ready to start.

**Treat with Veet® InfiniSilk™ Pro safely and properly!**  
Choose areas to treat that have not had recent exposure to the sun.

- The skin of the treatment areas should be shaved, clean, dry and free of any powders, antiperspirants or deodorants. Do not wax, pluck or tweeze the hair out. Shaving is an important step to get the results you want!
- Press the ON/STBY SWITCH. The POWER ON INDICATOR LIGHT (orange) will turn on and a fan

sound (similar to the sound of a hairdryer) will start.

- Approximately 3.5 seconds after pressing the ON/ STBY SWITCH, the READY INDICATOR LIGHT (green) will turn on. The device is then ready for you to trigger the first pulse at the lowest energy level.
- If this is your first treatment - Choose the energy level. InfiniSilk™ Pro offers 5 levels of energy from 1 (the lowest) to 5 (the highest). To choose the right energy level for your first treatment, test each treatment area separately:

- Place the TREATMENT SURFACE on the skin and press the trigger to deliver one pulse. If you did not experience abnormal discomfort press the + button once and deliver one pulse on a different spot at energy level 2.
- Continue testing at increased energy levels until you have reached the highest level that you feel comfortable with. If within one hour you experience no adverse effects, you can start your full treatment at that energy level.
- Repeat this test for each body part you intend to treat.

- Set the energy level. If the energy level should be higher than the lowest, press the ENERGY LEVEL SETTING BUTTONS using the "-" or "+" to respectively decrease or increase the energy level, until the desired energy level is set and represented by 3 energy level indicator lights (For details see "Energy Level" chapter inside this user manual).
- Apply the TREATMENT SURFACE to the skin, making sure the skin is spread evenly and smoothly. As soon as the TREATMENT SURFACE is in full contact with the skin, the READY INDICATOR LIGHT will start blinking.
- Press the PULSE BUTTON. The device will first determine the color of your skin. If the color of the skin is light enough for safe application, the device will flash a pulse of light onto your skin, and the READY INDICATOR LIGHT will turn off. You will see a bright flash of light and simultaneously hear a subtle pop sound, which is a normal noise for the device. You will feel a mild sensation of warmth and tingling. Veet® InfiniSilk™ Pro™ will immediately recharge for the next pulse. After 3.5 seconds the READY INDICATOR LIGHT will turn on again.
- Remove the TREATMENT SURFACE from the treated area of skin.
- Move the TREATMENT SURFACE to another area of skin. Use the pressure marks the device just made on your skin to guide you for proper positioning of the next pulse, avoiding both gaps and overlaps between pulses.

Veet® Infini'Silk™ Pro™ pulses should be administered in rows, starting at one end of each row and progressing sequentially towards the other end. This technique allows better control of skin coverage, and helps you avoid treating the same area more than once or overlapping skin areas. When applied to the skin, the Veet® Infini'Silk™ Pro™ TREATMENT SURFACE is designed to create temporary pressure marks on the treated area. These visible marks can be used for proper positioning of the next pulse.

Try to avoid overlapping pulses!

Warning: Do not treat the same area of skin more than once per hair removal session! Treating the same area of skin more than once per session increases the likelihood of adverse effects.

14. Repeat the process until the entire treatment area is covered.

If your skin blisters or burns, STOP USE IMMEDIATELY!

### Energy Level

Energy level determines the intensity of the Veet® Infini'Silk™ Pro™ light pulse delivered to your skin, from the lowest level (-) to the highest level (+). As energy level increases, so does hair removal results as well as the risk of possible side effects and complications.

Infini'Silk™ Pro offers 5 levels of energy from 1 (the lowest) to 5 (the highest), represented by 3 energy level indicator lights:

Level 1 (lowest) - the left indicator light is steadily lit

Level 2 - the left indicator light is steadily lit and the central indicator light is blinking.

Level 3 - both left and central indicator lights are steadily lit

Level 4 - the left and central indicator lights are steadily lit and the right indicator light is blinking.

Level 5 (highest) - all 3 indicator lights are steadily lit.

Whenever Veet® Infini'Silk™ Pro™ is turned on its energy level will automatically be reset to the

lowest energy level. Only one ENERGY LEVEL INDICATOR LIGHT will be on.

To set the energy level, press the ENERGY LEVEL SETTING BUTTONS using “-” or “+” to respectively decrease or increase the energy level.

### The Skin Color Sensor

Light-based hair removal on darker skin tones can result in adverse effects such as burns, blisters, and skin color changes (hyper- or hypo-pigmentation).

To prevent such misuse, a SKIN COLOR SENSOR in Veet® Infini'Silk™ Pro™ measures the treated skin complexion at the beginning of each session and occasionally during the session. If the SKIN COLOR SENSOR detects a skin tone that is too dark for Veet® Infini'Silk™ Pro™ application, the device will automatically stop emitting pulses.

If you see no light pulse and the SKIN TONE WARNING INDICATOR LIGHT blinks, this is an indication that your skin tone, measured by the SKIN COLOR SENSOR, is too dark for safe application. Try using the device on a different body part or contact RB Consumer Relations.

### Covering skin areas

Veet® Infini'Silk™ Pro™ pulses should be administered in rows, starting at one end of each row and progressing sequentially towards the other end. This technique allows better control of skin coverage, and helps you avoid treating the same area more than once or overlapping skin areas.

When applied to the skin, the Veet® Infini'Silk™ Pro™ TREATMENT SURFACE is designed to create temporary pressure marks on the treated area. These visible marks can be used for proper positioning of the next pulse.

Try to avoid overlapping pulses!

Do not treat the same area of skin more than once per hair removal session!

If your skin blisters or burns, STOP USE IMMEDIATELY!

### 9.1. What to Expect when treating with Veet® Infini'Silk™ Pro™?

For many people, using Veet® Infini'Silk™ Pro™ may be their first experience with a light-based device designed for home use. Veet® Infini'Silk™ Pro™ is simple to use, and hair removal sessions go by quickly. During a Veet® Infini'Silk™ Pro™ session it is normal to experience and feel:

- **A Fan Noise** – The cooling fan in Veet® Infini'Silk™ Pro™ makes noise similar to a hairdryer. This is normal.
- **A Pop Sound with Each Pulse** – When a pulse of HPL™ light is activated, it is normal to hear a subtle pop sound simultaneously with the flash of light.
- **Moderate Pressure of the Treatment Surface** – This is necessary and helpful for placement of adjacent pulses of light, and is part of the unique safety feature of Veet® Infini'Silk™ Pro™.
- **A Sensation of Warmth and Tingling** – During each pulse of light it is normal to feel a mild sensation of warmth and tingling from the light energy.
- **Some Mild Red or Pink Color** – During and just after your Veet® Infini'Silk™ Pro™ session it is not uncommon to see some very mild, pink-like color of the skin. This is usually most noticeable around the hairs themselves.

However if you see full redness of the skin, blistering or burns stop use of Veet® Infini'Silk™ Pro™ immediately.

### 9.2. After treating with Veet® Infini'Silk™ Pro™

- When Veet® Infini'Silk™ Pro™ session has been completed turn Veet® Infini'Silk™ Pro™ off by pressing the ON/STBY SWITCH. All the lit indicators will turn off and the STBY INDICATOR LIGHT will turn on (Be sure to remember the last energy level setting you used, as it will not be restored when turning Veet® Infini'Silk™ Pro™ on again.)
- Unplug the ADAPTOR cord from the electrical outlet.
- After each hair removal session it is recommended that you clean your Veet® Infini'Silk™ Pro device, especially the TREATMENT SURFACE (See: "Cleaning Veet® Infini'Silk™ Pro™").
- After cleaning, it is recommended to store your Veet® Infini'Silk™ Pro™ device in its original box, and keep it away from water.

#### Skin care following hair removal session

Do not expose treated areas of skin to the sun. Be sure to carefully protect the treated skin with sunscreen, throughout the hair removal period and for at least 2 weeks following the last Veet® Infini'Silk™ Pro™ session.

## 10. Maintenance of Veet® Infini'Silk™ Pro™

### 10.1. Cleaning Veet® Infini'Silk™ Pro

After each hair removal session, it is recommended to clean your Veet® Infini'Silk™ Pro™ device, and especially the TREATMENT SURFACE.

- Unplug Veet® Infini'Silk™ Pro™ before cleaning.
  - Use a dry, clean cloth and a specially formulated cleaner for electronic equipment to gently wipe Veet® Infini'Silk™ Pro™ surface, and especially the TREATMENT SURFACE.
- Never immerse Veet® Infini'Silk™ Pro™ or any of its parts in water!

### 10.2. Replacing the LAMP CARTRIDGE

The Veet® Infini'Silk™ Pro™ LAMP CARTRIDGE can fire a pre-set and limited number of light pulses. Pulse intensity is determined only according to the energy level setting of the device. There is no decrease of energy during the usable lifetime of the LAMP CARTRIDGE.

"When only 100 light pulses remain in the LAMP CARTRIDGE, the REPLACE CARTRIDGE INDICATOR LIGHT on the CONTROL PANEL will start blinking, indicating that the LAMP CARTRIDGE should soon be replaced."

"As the number of light pulses in the LAMP CARTRIDGE further decreases, the REPLACE CARTRIDGE INDICATOR LIGHT blinking will grow faster."

Once all the pulses in a LAMP CARTRIDGE have been used, the REPLACE CARTRIDGE INDICATOR LIGHT will remain steadily on, indicating that further pulses cannot be emitted.

In order to continue the hair removal session, it will now be necessary to replace the LAMP CARTRIDGE.

**Note:** A LAMP CARTRIDGE should also be replaced if large spots appear inside it or if its TREATMENT SURFACE is broken.

To replace a LAMP CARTRIDGE:

1. Press the ON/STBY SWITCH to turn Veet® Infini'Silk™ Pro™ to STBY mode.
2. Unplug the power adaptor from the electrical outlet.
3. Grasp the used LAMP CARTRIDGE on both sides, pull it out of the socket
4. Unwrap a new LAMP CARTRIDGE.
5. Push the new LAMP CARTRIDGE gently into place.

Note! After replacing the lamp cartridge with a new one, always set the energy to be one level lower than your normal setting. Raise the energy level back, only if you do not experience abnormal discomfort during or after the treatment.

## 11. Troubleshooting

### 11.1. “My Veet® InfiniSilk™ Pro™ does not start.”

- Make sure the POWER ADAPTOR is properly connected to the InfiniSilk™ Pro™ device.
- Make sure the POWER ADAPTOR is plugged into an electrical outlet on the wall.

### 11.2. “A light pulse is not emitted when I press the PULSE BUTTON

Make sure that the READY INDICATOR LIGHT is blinking to ensure you have good contact with the skin and that the TREATMENT SURFACE is evenly and firmly pressed to the skin. For your safety, the PULSE BUTTON will activate a pulse only if the TREATMENT SURFACE is firmly pressed against the skin.

- Check the SKIN TONE WARNING INDICATOR LIGHT. If it blinks, this is an indication that your skin tone, measured by the SKIN COLOR SENSOR, is too dark for safe application. Try using the device on a different body part or contact the RB® Consumer Relations
- Check the REPLACE CARTRIDGE INDICATOR LIGHT. If it is on, disconnect Veet® InfiniSilk™ Pro™ and replace the LAMP CARTRIDGE.
- Make sure that the READY INDICATOR is on.
  - If within 10 seconds the READY INDICATOR remains off turn Veet® InfiniSilk™ Pro™ off and back on by pressing the ON/STBY SWITCH twice.
  - If the problem persists, contact the RB Consumer Relations Center at 1-866-252-5327 or visit [www.veet.us](http://www.veet.us)
- Check the SYSTEM WARNING INDICATOR LIGHT.
  - If the light is on turn Veet® InfiniSilk™ Pro™ off and back on by pressing the ON/STBY SWITCH twice.
  - If the problem persists, contact the RB Consumer Relations Center at 1-866-252-5327 or visit [www.veet.us](http://www.veet.us)
- If the problem persists, contact the RB Consumer Relations Center at 1-866-252-5327 or visit [www.veet.us](http://www.veet.us)

Do not attempt to open or repair your Veet® InfiniSilk™ Pro device.

Opening Veet® InfiniSilk™ Pro™ may expose you to dangerous electrical components and to pulsed light energy, either of which may cause serious bodily damage and/or permanent eye injury. Trying to open Veet® InfiniSilk™ Pro™ may also damage the device and will void your warranty. Please contact RB Consumer Relations if you have a broken or damaged device in need or repair.

## 12. Consumer Relations

For more information about Veet® InfiniSilk™ or Veet products, contact the RB Consumer Relations Center at 1-866-252-5327 or visit [www.veet.us](http://www.veet.us)

If your device is broken, damaged, in need of repair, or for any other user assistance, contact the RB Consumer Relations Center at 1-866-252-5327 or visit [www.veet.us](http://www.veet.us).

## 13. Specifications

Spot Size	1.3cm x 3cm [3.9cm <sup>2</sup> ]	
Speed	1 pulse every 3.5 second: 1.1 cm <sup>2</sup> /sec	
Technology	Home Pulsed Light™	
Max Energy Level	Max 5J/cm <sup>2</sup>	
Wavelength	475-1200nm	
Charging Time/Power Source	Continuous operation	
Electrical Requirements	100-240VAC, 0.9A	
Operation and Safety	Skin Color Sensor ensures use only on suitable skin complexions.	
Package Size	22cm x 22cm x 14cm	
System weight	0.35kg	
Temperature	Operating	10 to 40°C
	Storage	-40 to +70°C
Relative Humidity	Operating	30 to 75%rH
	Storage	10 to 90%rH
Atmospheric pressure	Operating	700 to 1060hPa
	Storage	500 to 1060hPa

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